

Product datasheet for **TS403294P5**

PIAS2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PIAS2, transcript variant alpha, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC203294
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PIAS2 Mouse Monoclonal Antibody [Clone ID: OT15B4] (TA808334)
ACCN:	NM_173206 , NP_775298
Synonyms:	ARIP3; DIP; MIZ1; PIASX; SIZ2; ZMIZ4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_775298
Locus ID:	9063
Cytogenetics:	18q21.1
Protein Families:	Stem cell - Pluripotency, Stem cell relevant signaling - JAK/STAT signaling pathway, Transcription Factors



[View online »](#)

Protein Pathways: Jak-STAT signaling pathway, Pathways in cancer, Small cell lung cancer, Ubiquitin mediated proteolysis