

Product datasheet for **TS403346P5**

XTP4 (MIEN1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MIEN1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC203346
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	XTP4 (MIEN1) Mouse Monoclonal Antibody [Clone ID: OTI4D6] (TA504625)
ACCN:	NM_032339 , NP_115715
Synonyms:	C17orf37; C35; ORB3; RDX12; XTP4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_115715
Locus ID:	84299
Cytogenetics:	17q12



[View online »](#)