

## Product datasheet for **TS403397P5**

### Lin28 (LIN28A) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of LIN28A in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack   |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC203397  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)   |
| ACCN:                                 | <u>NM_024674</u> , <u>NP_078950</u>   |
| Synonyms:                             | CSDD1; LIN-28; lin-28A; LIN28; ZCCHC1   |
| Storage:                              | Room Temperature  |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <u>NP_078950</u>  |
| Locus ID:                             | 79727   |
| Cytogenetics:                         | 1p36.11   |
| Protein Families:                     | Transcription Factors   |



[View online »](#)