

Product datasheet for **TS403431P5**


PVRL4 (NECTIN4) CytoSection

Product data:


Product Type:	CytoSections
Description:	Transient overexpression of NECTIN4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC203431
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_030916</u> , <u>NP_112178</u>
Synonyms:	EDSS1; LNIR; nectin-4; PRR4; PVRL4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_112178</u>
Locus ID:	81607
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Adherens junction



[View online »](#)

Product images:

Immunohistology staining on NECTIN4 overexpressed cytosection [TS403431] with mouse anti DDK clone OTI11C3 C/N [TA180144] at 1:400 dilution 4C O/N. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC203431], 5 micron sections, 40x magnification



Immunohistology staining on PVRL4 overexpressed cytosection TS403431P5 with rabbit anti DDK clone OTIR5G2 C/N [TA592569] at 1:2000 dilution 5m RT. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone , 5 micron sections, 40x magnification