

Product datasheet for TS403506P5

ERK5 (MAPK7) CytoSection

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

| Product Type: | CytoSections |
|--|--|
| Description: | Transient overexpression of MAPK7, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC203506 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | ERK5 (MAPK7) Mouse Monoclonal Antibody [Clone ID: OTI5C7] (TA502203) |
| ACCN: | <u>NM 002749, NP 002740</u> |
| Synonyms: | BMK1; ERK4; ERK5; PRKM7 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP 002740</u> |
| Locus ID: | 5598 |
| Cytogenetics: | 17p11.2 |
| Protein Families: | Druggable Genome, Protein Kinase |
| Protein Pathways: | Gap junction, GnRH signaling pathway, MAPK signaling pathway, Neurotrophin signaling pathway |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US