

Product datasheet for **TS404072**

alpha Tubulin (TUBA3C) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TUBA3C, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204072
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_079836 , NP_524575
Synonyms:	bA408E5.3; TUBA2; tubulin, alpha 2; tubulin, alpha 3c
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_524575
Locus ID:	7278
Cytogenetics:	13q12.11
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection



[View online »](#)