

Product datasheet for **TS404396**

GABPB2 (GABPB1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GABPB1, transcript variant gamma-1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204396
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GABPB2 (GABPB1) Mouse Monoclonal Antibody [Clone ID: OTI7A7] (TA811268)
ACCN:	<u>NM_002041</u> , <u>NP_002032</u>
Synonyms:	BABPB2; E4TF1; E4TF1-47; E4TF1-53; E4TF1B; GABPB; GABPB-1; GABPB2; NRF2B1; NRF2B2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_002032</u>
Locus ID:	2553
Cytogenetics:	15q21.2
Protein Families:	Transcription Factors



[View online »](#)