

Product datasheet for **TS404535**

Neurotrophin 4 (NTF4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NTF4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204535
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Neurotrophin 4 (NTF4) Mouse Monoclonal Antibody [Clone ID: OTI2E8] (TA500034)
ACCN:	NM_006179 , NP_006170
Synonyms:	GLC10; GLC10; NT-4; NT-4/5; NT-5; NT4; NT5; NTF5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006170
Locus ID:	4909
Cytogenetics:	19q13.33
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	MAPK signaling pathway, Neurotrophin signaling pathway



[View online »](#)