

Product datasheet for **TS404546P5**

GEF H1 (ARHGEF2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ARHGEF2, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204546
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GEF H1 (ARHGEF2) Mouse Monoclonal Antibody [Clone ID: OTI5F2] (TA809320)
ACCN:	NM_004723 , NP_004714
Synonyms:	GEF; GEF-H1; GEFH1; Lfc; LFP40; NEDMHM; P40
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004714
Locus ID:	9181
Cytogenetics:	1q22
Protein Families:	Druggable Genome
Protein Pathways:	Pathogenic Escherichia coli infection



[View online »](#)