

## Product datasheet for **TS404711P5**

### Caspase 4 (CASP4) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CASP4, transcript variant gamma, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204711
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Caspase 4 (CASP4) Mouse Monoclonal Antibody [Clone ID: OTI4A2] (TA805590)
ACCN:	<a href="#">NM_033306</a> , <a href="#">NP_150649</a>
Synonyms:	ICE(rel)II; ICEREL-II; ICH-2; Mih1; Mih1/TX; TX
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_150649</a>
Locus ID:	837
Cytogenetics:	11q22.3
Protein Families:	Druggable Genome, Protease



[View online »](#)