

Product datasheet for **TS404763**

GNPDA1 CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of GNPDA1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC204763 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | NM_005471 , NP_005462 |
| Synonyms: | GNP1; GNPDA; GNPI; GPI; HLN |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_005462 |
| Locus ID: | 10007 |
| Cytogenetics: | 5q31.3 |
| Protein Pathways: | Amino sugar and nucleotide sugar metabolism, Metabolic pathways |



[View online »](#)