

Product datasheet for **TS404895P5**

p15 INK4b (CDKN2B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CDKN2B, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204895
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	p15 INK4b (CDKN2B) Mouse Monoclonal Antibody [Clone ID: OTI3B6] (TA505206)
ACCN:	NM_004936 , NP_004927
Synonyms:	CDK4I; INK4B; MTS2; P15; p15INK4b; TP15
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004927
Locus ID:	1030
Cytogenetics:	9p21.3
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Pathways in cancer, Small cell lung cancer, TGF-beta signaling pathway



[View online »](#)