

Product datasheet for **TS404904**

SERPINE2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SERPINE2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204904
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SERPINE2 Mouse Monoclonal Antibody [Clone ID: OT12C9] (TA504611)
ACCN:	NM_006216 , NP_006207
Synonyms:	GDN; GDNPF; PI-7; PI7; PN-1; PN1; PNI
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006207
Locus ID:	5270
Cytogenetics:	2q36.1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)