

## Product datasheet for **TS404924P5**

### alpha Actinin (ACTN1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACTN1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC204924
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	alpha Actinin (ACTN1) Mouse Monoclonal Antibody [Clone ID: OT17A4] (TA500072)
ACCN:	<u><a href="#">NM_001102</a></u> , <u><a href="#">NP_001093</a></u>
Synonyms:	BDPLT15
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001093</a></u>
Locus ID:	87
Cytogenetics:	14q22-q24
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS



[View online »](#)

**Protein Pathways:**

Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Focal adhesion, Leukocyte transendothelial migration, Regulation of actin cytoskeleton, Systemic lupus erythematosus, Tight junction