

Product datasheet for **TS405035P5**

ZKSCAN1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZKSCAN1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205035
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ZKSCAN1 Mouse Monoclonal Antibody [Clone ID: OTI2C9] (TA810966)
ACCN:	<u>NM_003439</u> , <u>NP_003430</u>
Synonyms:	KOX18; PHZ-37; ZNF36; ZNF139; ZSCAN33
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003430</u>
Locus ID:	7586
Cytogenetics:	7q22.1
Protein Families:	Transcription Factors


[View online »](#)