

Product datasheet for **TS405248P5**

TM4SF2 (TSPAN7) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TSPAN7 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205248
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	TM4SF2 (TSPAN7) Mouse Monoclonal Antibody [Clone ID: OTI2D6] (TA811638)
ACCN:	NM_004615 , NP_004606
Synonyms:	A15; CCG-B7; CD231; DXS1692E; MRX58; MXS1; TALLA-1; TM4SF2; TM4SF2b
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004606
Locus ID:	7102
Cytogenetics:	Xp11.4
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)