

## Product datasheet for **TS405289**

### G protein alpha inhibitor 1 (GNAI1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GNAI1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205289
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	G protein alpha inhibitor 1 (GNAI1) Mouse Monoclonal Antibody [Clone ID: OTI2D1] (TA811907)
ACCN:	<u><a href="#">NM_002069</a></u> , <u><a href="#">NP_002060</a></u>
Synonyms:	Gi
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_002060</a></u>
Locus ID:	2770
Cytogenetics:	7q21.11
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:**

Axon guidance, Chemokine signaling pathway, Gap junction, Leukocyte transendothelial migration, Long-term depression, Melanogenesis, Progesterone-mediated oocyte maturation, Tight junction