

Product datasheet for **TS405355P5**

USP12 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of USP12 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205355
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	USP12 Mouse Monoclonal Antibody [Clone ID: OTI2A6] (TA800037)
ACCN:	NM_182488 , NP_872294
Synonyms:	UBH1; USP12L1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_872294
Locus ID:	219333
Cytogenetics:	13q12.13
Protein Families:	Druggable Genome, Protease



[View online »](#)