

Product datasheet for **TS405372P5**

PYK2 (PTK2B) CytoSection

Product data:

| | |
|--|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of PTK2B, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC205372 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | <u>NM_173174</u> , <u>NP_775266</u> |
| Synonyms: | CADTK; CAKB; FADK2; FAK2; PKB; PTK; PYK2; RAFTK |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_775266</u> |
| Locus ID: | 2185 |
| Cytogenetics: | 8p21.2 |
| Protein Families: | Druggable Genome, Protein Kinase |
| Protein Pathways: | Calcium signaling pathway, Chemokine signaling pathway, GnRH signaling pathway, Leukocyte transendothelial migration, Natural killer cell mediated cytotoxicity |



[View online »](#)