

Product datasheet for **TS405447P5**

ATP6V1B2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ATP6V1B2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205447
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	ATP6V1B2 Mouse Monoclonal Antibody [Clone ID: OTI1E11] (TA505244)
ACCN:	NM_001693 , NP_001684
Synonyms:	ATP6B1B2; ATP6B2; DOOD; HO57; VATB; Vma2; VPP3; ZLS2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001684
Locus ID:	526
Cytogenetics:	8p21.3
Protein Families:	Druggable Genome
Protein Pathways:	Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection



[View online »](#)