

Product datasheet for **TS405621P5**

UBE2D1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UBE2D1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205621
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UBE2D1 Mouse Monoclonal Antibody [Clone ID: OT11E7] (TA803633)
ACCN:	NM_003338 , NP_003329
Synonyms:	E2(17)KB1; SFT; UBC4/5; UBCH5; UBCH5A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003329
Locus ID:	7321
Cytogenetics:	10q21.1
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)