

Product datasheet for **TS405843**

FARS2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FARS2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205843
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_006567 , NP_006558
Synonyms:	COXPD14; FARS1; HSPC320; mtPheRS; PheRS; SPG77
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006558
Locus ID:	10667
Cytogenetics:	6p25.1
Protein Pathways:	Aminoacyl-tRNA biosynthesis



[View online »](#)