

## Product datasheet for **TS405916P5**

### ST2 (IL1RL1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL1RL1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205916
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_003856</a> , <a href="#">NP_003847</a>
Synonyms:	DER4; FIT-1; IL33R; ST2; ST2L; ST2V; T1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003847</a>
Locus ID:	9173
Cytogenetics:	2q12.1
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)