

Product datasheet for **TS405920P5**

Neurofilament (NEFL) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NEFL in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205920
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Neurofilament (NEFL) Mouse Monoclonal Antibody [Clone ID: OTI11F6] (TA801192)
ACCN:	NM_006158 , NP_006149
Synonyms:	CMT1F; CMT2E; CMTDIG; NF-L; NF68; NFL; PPP1R110
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006149
Locus ID:	4747
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS
Protein Pathways:	Amyotrophic lateral sclerosis (ALS)



[View online »](#)