

Product datasheet for **TS405940P5**

ERp57 (PDIA3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PDIA3 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205940
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	ERp57 (PDIA3) Mouse Monoclonal Antibody [Clone ID: OTI3H5] (TA505010)
ACCN:	NM_005313 , NP_005304
Synonyms:	ER60; ERp57; ERp60; ERp61; GRP57; GRP58; HEL-S-93n; HEL-S-269; HsT17083; P58; PI-PLC
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005304
Locus ID:	2923
Cytogenetics:	15q15.3
Protein Families:	Druggable Genome
Protein Pathways:	Antigen processing and presentation



[View online »](#)