

Product datasheet for **TS405943P5**

Albumin (ALB) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ALB in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC205943
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Albumin (ALB) Mouse Monoclonal Antibody [Clone ID: OTI3B10] (TA507286)
ACCN:	NM_000477 , NP_000468
Synonyms:	FDAHT; HSA; PRO0883; PRO0903; PRO1341
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000468
Locus ID:	213
Cytogenetics:	4q13.3
Protein Families:	Secreted Protein



[View online »](#)