

## Product datasheet for **TS406173**

### RAB11FIP5 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RAB11FIP5 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206173
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_015470</a></u> , <u><a href="#">NP_056285</a></u>
Synonyms:	<u><a href="#">gaf-1</a></u> ; <u><a href="#">GAF1</a></u> ; <u><a href="#">pp75</a></u> ; <u><a href="#">rab11-FIP5</a></u> ; <u><a href="#">RIP11</a></u>
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_056285</a></u>
Locus ID:	26056
Cytogenetics:	2p13.2
Protein Pathways:	Endocytosis



[View online »](#)