

Product datasheet for **TS406335P5**

U2AF1L3 (U2AF1L4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of U2AF1L4, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206335
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	U2AF1L3 (U2AF1L4) Mouse Monoclonal Antibody [Clone ID: OTI8F10] (TA810408)
ACCN:	NM_144987 , NP_659424
Synonyms:	U2AF1-RS3; U2AF1L3; U2AF1L3V1; U2AF1RS3; U2af26
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_659424
Locus ID:	199746
Cytogenetics:	19q13.12



[View online »](#)