

Product datasheet for **TS406457P5**

Rad51L1 (RAD51B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RAD51B, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206457
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Rad51L1 (RAD51B) Mouse Monoclonal Antibody [Clone ID: OTI7H4] (TA503550)
ACCN:	NM_133509 , NP_598193
Synonyms:	R51H2; RAD51L1; REC2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_598193
Locus ID:	5890
Cytogenetics:	14q24.1
Protein Families:	Druggable Genome
Protein Pathways:	Homologous recombination



[View online »](#)