

Product datasheet for **TS406535P5**

Estrogen Sulfotransferase (SULT1E1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SULT1E1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206535
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_005420 , NP_005411
Synonyms:	EST; EST-1; ST1E1; STE
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005411
Locus ID:	6783
Cytogenetics:	4q13.3
Protein Pathways:	Androgen and estrogen metabolism, Sulfur metabolism



[View online »](#)