

## Product datasheet for **TS406640P5**

### Pyruvate Dehydrogenase E2 (DLAT) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DLAT in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206640
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_001931</a> , <a href="#">NP_001922</a>
Synonyms:	DLTA; E2; PBC; PDC-E2; PDCE2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001922</a>
Locus ID:	1737
Cytogenetics:	11q23.1
Protein Families:	Druggable Genome
Protein Pathways:	Citrate cycle (TCA cycle), Glycolysis / Gluconeogenesis, Metabolic pathways, Pyruvate metabolism



[View online »](#)