

## Product datasheet for **TS406793**

### ABAT CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ABAT, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206793
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ABAT Mouse Monoclonal Antibody [Clone ID: OTI3D2] (TA806896)
ACCN:	<a href="#">NM_000663</a> , <a href="#">NP_000654</a>
Synonyms:	GABA-AT; GABAT; NPD009
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_000654</a>
Locus ID:	18
Cytogenetics:	16p13.2
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:**

Alanine, aspartate and glutamate metabolism, beta-Alanine metabolism, Butanoate metabolism, Metabolic pathways, Propanoate metabolism, Valine, leucine and isoleucine degradation