

Product datasheet for **TS406830**

DHLAG (CD74) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CD74, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC206830
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DHLAG (CD74) Mouse Monoclonal Antibody [Clone ID: OTI1H3] (TA507339)
ACCN:	NM_001025158 , NP_001020329
Synonyms:	DHLAG; HLADG; Ia-GAMMA; II; p33
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001020329
Locus ID:	972
Cytogenetics:	5q33.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Antigen processing and presentation



[View online »](#)