

Product datasheet for TS406855P5

ACADS CytoSection

Product data:

Product Type: CytoSections **Description:** Transient overexpression of ACADS in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack Species: Human **Expression Host:** HEK293T **Expression cDNA Clone** TrueORF Clone RC206855 or AA Sequence: C-MYC/DDK Tag: **Detection Antibodies:** Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) **Target Detection** ACADS Mouse Monoclonal Antibody [Clone ID: OTI1D2] (TA800429) Antibodies: ACCN: NM 000017, NP 000008 ACAD3; SCAD Synonyms: Storage: **Room Temperature** Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. **Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. This product is for research use only and is not approved for use in humans or in clinical Note: diagnosis. **RefSeq:** NP 000008 35 Locus ID: **Cytogenetics:** 12q24.31 **Protein Families:** Druggable Genome **Protein Pathways:** Butanoate metabolism, Fatty acid metabolism, Metabolic pathways, Valine, leucine and isoleucine degradation



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn