

## Product datasheet for **TS407020**

### IL10RA CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL10RA, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207020
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IL10RA Mouse Monoclonal Antibody [Clone ID: OT11D10] (TA506929)
ACCN:	<a href="#">NM_001558</a> , <a href="#">NP_001549</a>
Synonyms:	CD210; CD210a; CDW210A; HIL-10R; IL-10R1; IL10R
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001549</a>
Locus ID:	3587
Cytogenetics:	11q23.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway



[View online »](#)