

## Product datasheet for **TS407102P5**

### Hsp105 (HSPH1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HSPH1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207102
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_006644</a></u> , <u><a href="#">NP_006635</a></u>
Synonyms:	HSP105; HSP105A; HSP105B; NY-CO-25
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_006635</a></u>
Locus ID:	10808
Cytogenetics:	13q12.3
Protein Families:	Stem cell - Pluripotency



[View online »](#)