

Product datasheet for **TS407107P5**

DHR SX CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DHR SX in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207107
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_145177</u> , <u>NP_660160</u>
Synonyms:	CXorf11; DHR S5X; DHR S5Y; DHR SX Y; DHR SY; SDR7C6; SDR46C1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_660160</u>
Locus ID:	207063
Cytogenetics:	X;Y
Protein Families:	Druggable Genome



[View online »](#)