

## Product datasheet for **TS407138P5**

### Lipin 1 (LPIN1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LPIN1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207138
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Lipin 1 (LPIN1) Mouse Monoclonal Antibody [Clone ID: OTI8F10] (TA805972)
ACCN:	<u><a href="#">NM_145693</a></u> , <u><a href="#">NP_663731</a></u>
Synonyms:	PAP1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_663731</a></u>
Locus ID:	23175
Cytogenetics:	2p25.1



[View online »](#)