

Product datasheet for **TS407276P5**

Monoamine Oxidase A (MAOA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAOA in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207276
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Monoamine Oxidase A (MAOA) Mouse Monoclonal Antibody [Clone ID: OT12F10] (TA504266)
ACCN:	NM_000240 , NP_000231
Synonyms:	BRNRS; MAO-A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000231
Locus ID:	4128
Cytogenetics:	Xp11.3
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways:

Arginine and proline metabolism, Drug metabolism - cytochrome P450, Glycine, serine and threonine metabolism, Histidine metabolism, Metabolic pathways, Phenylalanine metabolism, Tryptophan metabolism, Tyrosine metabolism