

Product datasheet for **TS407384P5**

Kv beta 1 (KCNAB1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNAB1, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207384
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Kv beta 1 (KCNAB1) Mouse Monoclonal Antibody [Clone ID: OTI7F12] (TA503971)
ACCN:	NM_172159 , NP_751891
Synonyms:	AKR6A3; hKvb3; hKvBeta3; KCNA1B; KV-BETA-1; Kvb1.3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_751891
Locus ID:	7881
Cytogenetics:	3q25.31
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)