

## Product datasheet for **TS407414P5**

### MAPK4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAPK4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207414
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MAPK4 Mouse Monoclonal Antibody [Clone ID: OT12F4] (TA505872)
ACCN:	<a href="#">NM_002747</a> , <a href="#">NP_002738</a>
Synonyms:	ERK-4; ERK4; p63-MAPK; p63MAPK; PRKM4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_002738</a>
Locus ID:	5596
Cytogenetics:	18q21.1-q21.2
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)