

## Product datasheet for **TS407522P5**

### ADAR1 (ADAR) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of ADAR, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC207522  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| ACCN:                                 | <u><a href="#">NM_001111</a></u> , <u><a href="#">NP_001102</a></u>   |
| Synonyms:                             | ADAR1; AGS6; DRADA; DSH; DSRAD; G1P1; IFI-4; IFI4; K88DSRBP; P136   |
| Storage:                              | Room Temperature  |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <u><a href="#">NP_001102</a></u>  |
| Locus ID:                             | 103   |
| Cytogenetics:                         | 1q21.3  |
| Protein Families:                     | Druggable Genome  |
| Protein Pathways:                     | Cytosolic DNA-sensing pathway   |



[View online »](#)