

Product datasheet for **TS407568P5**

RNF41 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RNF41, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207568
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RNF41 Mouse Monoclonal Antibody [Clone ID: OTI2C1] (TA811158)
ACCN:	NM_194358 , NP_919339
Synonyms:	FLRF; NRDP1; SBB103
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_919339
Locus ID:	10193
Cytogenetics:	12q13.3
Protein Pathways:	Endocytosis



[View online »](#)