

Product datasheet for **TS407577P5**

Steroidogenic Factor 1 (NR5A1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NR5A1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207577
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Steroidogenic Factor 1 (NR5A1) Mouse Monoclonal Antibody [Clone ID: OT11B4] (TA806157)
ACCN:	NM_004959 , NP_004950
Synonyms:	AD4BP; ELP; FTZ1; FTZF1; hSF-1; POF7; SF-1; SF1; SPGF8; SRXX4; SRXY3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004950
Locus ID:	2516
Cytogenetics:	9q33.3
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)