

## Product datasheet for **TS407596**

### TRAIL (TNFSF10) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TNFSF10, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207596
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TRAIL (TNFSF10) Mouse Monoclonal Antibody [Clone ID: OT112E12] (TA811771)
ACCN:	<a href="#">NM_003810</a> , <a href="#">NP_003801</a>
Synonyms:	Apo-2L; APO2L; CD253; TL2; TNLG6A; TRAIL
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003801</a>
Locus ID:	8743
Cytogenetics:	3q26.31
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Apoptosis, Cytokine-cytokine receptor interaction, Natural killer cell mediated cytotoxicity



[View online »](#)