

## Product datasheet for **TS407599P5**

### TFF1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TFF1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207599
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TFF1 Mouse Monoclonal Antibody [Clone ID: OT11A4] (TA804518)
ACCN:	<a href="#">NM_003225</a> , <a href="#">NP_003216</a>
Synonyms:	BCEI; D21S21; HP1.A; HPS2; pNR-2; pS2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_003216</a>
Locus ID:	7031
Cytogenetics:	21q22.3
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)