

Product datasheet for **TS407616**

NFkB Inducing Kinase NIK (MAP3K14) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAP3K14 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207616
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_003954 , NP_003945
Synonyms:	FTDCR1B; HS; HSNIK; NIK
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003945
Locus ID:	9020
Cytogenetics:	17q21.31
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Apoptosis, Epithelial cell signaling in Helicobacter pylori infection, MAPK signaling pathway, T cell receptor signaling pathway



[View online »](#)