

Product datasheet for **TS407638P5**

DOLPP1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DOLPP1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207638
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_020438</u> , <u>NP_065171</u>
Synonyms:	LSFR2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_065171</u>
Locus ID:	57171
Cytogenetics:	9q34.11
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	N-Glycan biosynthesis



[View online »](#)