

Product datasheet for **TS407761P5**

ST3GAL2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ST3GAL2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207761
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_006927 , NP_008858
Synonyms:	Gal-NAc6S; SIAT4B; ST3GalA.2; ST3GALII
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_008858
Locus ID:	6483
Cytogenetics:	16q22.1
Protein Families:	Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - ganglio series, Glycosphingolipid biosynthesis - globo series, Keratan sulfate biosynthesis, Metabolic pathways, O-Glycan biosynthesis



[View online »](#)