

Product datasheet for **TS407780**

Peroxiredoxin 6 (PRDX6) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PRDX6 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC207780
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Peroxiredoxin 6 (PRDX6) Mouse Monoclonal Antibody [Clone ID: OT14D1] (TA813653)
ACCN:	NM_004905 , NP_004896
Synonyms:	1-Cys; aiPLA2; AOP2; HEL-S-128m; LPCAT-5; NSGPx; p29; PRX
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004896
Locus ID:	9588
Cytogenetics:	1q25.1
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Methane metabolism, Phenylalanine metabolism



[View online »](#)